

Summary of Zones of Entrapment Dimensions as Recognized by the F.D.A. verses SleepSafe Beds' Dimensioning

EXHIBIT ONE:

EXHIBIT ONE:		FDA Standards	New York State Standards	SleepSafe™ Actuals
Zone 1	Within the Rail	4 3/4"	2 1/3"	Zero"
Zone 2	Between the Top of the Compressed Mattress and the Bottom of the Rail Between Rail Supports: see below *	4 3/4"	Zero"	Zero " See Zone 3 and "Zone B" NO Lateral Play in Guard Rail.
Zone 3	Between the Rail and the Mattress	4 3/4"	Open hand Approx. 3/4"	Width of Open Hand, Approx 3/4" NO Lateral Play in Guard Rail!
Zone 4	Between the Top of the Compressed Mattress and the Bottom of the Rail at the End of the Rail: see below *	2 1/3"	Does not recognize this Zone	Zero " See Zone 3 and "Zone B" NO Lateral Play in Guard Rail.
Zone 5	Between the Split Bed Rails	2 1/3" -- 12 1/2"	2 13" -- 10"	Zero" / Not Applicable
Zone 6	Between the End of the Rail and the Side Edge of the Head or Foot Board	Head: 2 1/3" Foot: 12."	Head: 2 1/3" Foot: 10"	Head: 3/4" Foot: 3/4"
Zone 7	Between the Head or Foot Board and the End of the Mattress	4 3/4"	Open Hand Approx. 3/4"	Zero" or Tight Fit

DRAFT GUIDANCE for INDUSTRY and FDA STAFF:

"Zone A"	Between Sleeping Surface and Bottom of Head or Foot Board	<u>Not Specified</u>	NS	> 8" Overlap
"Zone B"	* Between FLAT Sleeping Surface and Bottom of Guard Rail	NS	NS	4" Overlap

ALL THE ABOVE ARE IDENTIFIED WITH DOCKET NUMBER 2004-0343.